

## Vidre



# Vidoire



Prohs Vidoire purpose is the ousting of organic substances in liquid or solid state (ex.: organic fluids such as escarpments or urine), eliminating them in a hygienic and secure form, guaranteeing its proper processing.

## Technical Characteristics

- Fully built in AISI 304 stainless steel.
- Mobile grating for containers support, to allow a more secure usage.
- Siphon in silicone, with wide passage section (diam. 90 mm), with discharge to the ground or wall.
- Connection to public water supply.



Brand	Model	Dimensions (mm) W. x L. x H.
PROHS	P	470 x 430 x 595
	G	470 x 560 x 595

The flusher must be placed at 1100mm from the ground (above the vidoire).

## Option

- With cover

Due to product improvement and development, PROHS reserves the right to make changes without notice.