

**PROHS** Equipamento Hospitalar  
e Serviços Associados, S.A.

**Rua do Castanhal, 316**  
Z. Industrial Maia I Sector II  
Apartado 6019 EC Outeiro  
4476-908 **Maia — Portugal**

| **T** | +351 229 059 170  
| **F** | +351 229 015 900  
| **E** | [prohs@prohs.pt](mailto:prohs@prohs.pt)  
| **W** | [www.prohs.pt](http://www.prohs.pt)

# Wall and Central Rack



**PROHS®**

# Wall and Central Rack

Equipment that is intended for the warehouse of Sterilization. They allow the proper storage of sterilized material.



## Wall Rack

### Technical Characteristics

- Fully built in AISI 304 stainless steel.
- Special system designed for wired baskets.
- Allows the storage of baskets with several sizes.
- Dimensions (height x thickness): 2000 x 50mm

**Can be produced in different dimensions.**



## Central Rack

### Technical Characteristics

- Fully built in AISI 304 stainless steel.
- They have a support structure itself, allowing to storage baskets on both sides.
- Special system designed for wired baskets.
- Allows the storage of baskets with various measures.
- Dimensions (width x height): 1700 x 2000 mm (3 sections model)

**Can be produced in different dimensions.**

### Option

- Revolving casters with or without brake.

Due to product improvement and development, PROHS reserves the right to make changes without notice.