

PROHS Equipamento Hospitalar
e Serviços Associados, S.A.

Rua do Castanhal, 316
Z. Industrial Maia I Sector II
Apartado 6019 EC Outeiro
4476-908 **Maia — Portugal**

| **T** | +351 229 059 170
| **F** | +351 229 015 900
| **E** | prohs@prohs.pt
| **W** | www.prohs.pt

Shelved Stand

PROHS[®]



Shelved Stand

Meant for the storing of sterilized material, accommodated in trays, boxes, containers, and others.
Can be used in various areas of an hospital.



Stainless Steel Stand with shelves

Technical Characteristics

- Fully built in AISI 304L stainless steel.
- With 5 shelves.
- Brushed finish.
- Lightweight.
- Easy storage and access.
- Regulable feet, for every surface.
- Dimensions (WxLxH): 900x400x2000mm

Can be produced in different dimensions.

Option

- Height adjustable shelves.



Stainless Steel Stand with wired shelves

Technical Characteristics

- Fully built in AISI 304L stainless steel.
- With 5 shelves.
- Glazed finish.
- Lightweight.
- Easy storage and access.
- Regulable feet, for every surface.
- Dimensions (WxLxH): 900x400x2000mm

Can be produced in different dimensions.

Option

- Four rotating wheels (two with brakes).

Due to product improvement and development, PROHS reserves the right to make changes without notice.