

**PROHS** Equipamento Hospitalar  
e Serviços Associados, S.A.

**Rua do Castanhal, 316**  
Z. Industrial Maia I Sector II  
Apartado 6019 EC Outeiro  
4476-908 **Maia — Portugal**

**| T |** +351 229 059 170  
**| F |** +351 229 015 900  
**| E |** [prohs@prohs.pt](mailto:prohs@prohs.pt)  
**| W |** [www.prohs.pt](http://www.prohs.pt)

# Cabinet / Bedpan Support

**PROHS®**



# Cabinet / Bedpan Support

Cabinet / Wall Support for storage of bedpans and urinals after its disinfection.




## Cabinet

### Technical Characteristics

- Built in AISI 304 stainless steel.
- Capacity to 6 bedpans
- Dimensions: 800 x 400 x 620 mm (W\*D\*H)

**Can be manufactured in different sizes.**



## Suporte

### Technical Characteristics

- Built in AISI 304 stainless steel.
- Capacity to 6 bedpans and 6 urinals.
- Dimensions: 800 x 300 x 530 mm (W\*D\*H)

**Can be manufactured in different sizes.**

Due to product improvement and development, PROHS reserves the right to make changes without notice.