

Preparation and Packing Table

PROHS Equipamento Hospitalar
e Serviços Associados, S.A.

Rua do Castanhal, 316
Z. Industrial Maia I Sector II
Apartado 6019 EC Outeiro
4476-908 Maia — Portugal

| T | +351 229 059 170
| F | +351 229 015 900
| E | prohs@prohs.pt
| W | www.prohs.pt



Preparation and Packing Table

The Preparation and Packing Tables are specially developed for the sterilization rooms.

Technical Characteristics

- Made in stainless steel AISI 304L with flat tops and soundproofed, the support consists of a welded tubular Argon structures.
- Rounded and polished finish for easy cleaning and promoting hygiene products.
- Structure of square tubing in stainless steel AISI 304, with adjustable feet height.
- The top shelf, tilted for easy access by the operator. Shelf less than the length of the table.
- Accommodate various accessories: paper holder, holder baskets, drawers, etc.
- Possibility of having two opposite positions.

Can be produced in different dimensions.



Preparation and Packing Table - double post
2000 x 850 x 900mm



Preparation and Packing Table - simple
2000 x 1000 x 900mm



Preparation and Packing Table - simple post
1500 x 700 x 900mm

Due to product improvement and development, PROHS reserves the right to make changes without notice.